

Laboratory Specimen No: 70505981
Client Account: 10932 / RM001Q
Location / Contact: HARRISON
Report To: Rodger's Memorial Hospital ATTN: Jim Smith, DO 456 Main Street Anytown, IN 46001 Fx: 812-555-1213

Donor's Name:	DOE, JOHN	
Donor ID/SSN:	123-89-4567	
Specimen ID Number:	RR1654978	
D.O.B. / Gender:	04/23/1978	Male
Reason for Testing:	Not given	
Date / Time Collected:	05/10/2010	13:50:00
Date Received:	05/11/2010	
Date Reported:	05/13/2010	

Container(s): 01: Urine Urine, Random	Test(s): 70195	DetectiMed Panel, Expanded (Urine)
---	-----------------------	------------------------------------

Analyte Name	Result	Cut-off	Unit	Therapeutic	Loc
AMPHETAMINES	Negative	500	ng/mL		
BARBITURATES	Negative	200	ng/mL		
BENZODIAZEPINES	POSITIVE	75	ng/mL		
Lorazepam	POSITIVE	50	ng/mL		
Lorazepam, Quant		126	ng/mL		
BUPRENORPHINE/METABOLITE	Negative	5	ng/mL		
CANNABINOIDS	Negative	20	ng/mL		
CARISOPRODOL/METABOLITE	Negative	100	ng/mL		
COCAINE/METABOLITES	POSITIVE	150	ng/mL		
Benzoyllecgonine	POSITIVE	30	ng/mL		
Benzoyllecgonine, Quant		48	ng/mL		
FENTANYL	Negative	2.0	ng/mL		
METHADONE/METABOLITE	POSITIVE	100	ng/mL		
Methadone	POSITIVE	100	ng/mL		
Methadone, Quant		5000	ng/mL		
EDDP	NEGATIVE	100	ng/mL		
OPIATES	Negative	50	ng/mL		
OXYCODONE/METABOLITE	Negative	50	ng/mL		
PROPOXYPHENE/METABOLITE	Negative	300	ng/mL		
ALCOHOLS	Negative	0.02	% (w/v)		
TRAMADOL/METABOLITE	Negative	200	ng/mL		
ACETAMINOPHEN	Negative	10	ug/mL		
CREATININE	Normal	5	mg/dL		
Creatinine, Quant		35	mg/dL		
pH		6.7		4.50-9.00	
GENERAL OXIDANTS	Negative	200	ug/mL		
MEDICATIONS LISTED:	=====				
1)	Methadone				
2)	Lorazepam				

Laboratory Specimen No: 70505981

Donor's Name: DOE, JOHN

Urine Drug Test – Quick Reference

“Positive Consistent” drug (or metabolite) detected and IS included in list of medications
 “Positive Inconsistent” drug (or metabolite) detected and is NOT included in list of medications
 “Negative Inconsistent” drug (or metabolite) NOT detected or below cutoff level and IS included in list of medications

Benzodiazepines

Drug	Reason for Positive Consistent / Inconsistent	Detection Window
Lorazepam	Use of Lorazepam (Ativan)	1-10 Days

Cocaine

Drug	Reason for Positive Consistent / Inconsistent	Detection Window
Benzoylcegonine	Metabolite of Cocaine	1-4 Days

Methadone

Drug	Reason for Positive Consistent / Inconsistent	Detection Window
Methadone	Use of Methadone (Dolophine, Methadose)	1-10 Days (pH dependent)

DISCLAIMER: This is a reference document that outlines the content and formatting traditionally included on AIT client reports. For current information about testing capabilities, panel options, and cut-off levels, please contact Client Services at 800 875-3894.

If you have any additional questions, feel free to contact us at: toxicologist@aitlabs.com